


# Turan<sup>TM</sup>

## URETERAL ACCESS SHEATH




- PTFE liner facilitates the entry of ureteroscopes and instruments
- Kink Resistance: Sheath is reinforced with stainless steel coil, designed provide stability and strength during placement
- Radiopaque stainless coil facilitates placement of the sheath under X-Ray
- Hydrophilic coated sheath offers lubricity during placement minimizing the risk of trauma

Turan™ ureteral access sheath provides dilation and working channel for the entry of ureteroscopes and instrument. It is intended to protect the ureter from trauma during repeated instrument exchanges during ureteroscopy procedures.

## Specification

Order No.	Size (ID Fr*Length )	Effective Length	Internal Diameter	External Diameter
		(mm)±5%	(mm)±0.15	(mm)±0.15
71026501	10Fr*265mm	265	3.3	4
71036501	10Fr*365mm	365	3.3	4
71041501	10Fr*415mm	415	3.3	4
71046501	10Fr*465mm	465	3.3	4
71056501	10Fr*565mm	565	3.3	4
71226501	12Fr*265mm	265	4.1	4.7
71236501	12Fr*365mm	365	4.1	4.7
71241501	12Fr*415mm	415	4.1	4.7
71246501	12Fr*465mm	465	4.1	4.7
71256501	12Fr*565mm	565	4.1	4.7
71326501	13Fr*265mm	265	4.35	5
71336501	13Fr*365mm	365	4.35	5
71341501	13Fr*415mm	415	4.35	5
71346501	13Fr*465mm	465	4.35	5
71356501	13Fr*565mm	565	4.35	5
71426501	14Fr*265mm	265	4.7	5.3
71436501	14Fr*365mm	365	4.7	5.3
71441501	14Fr*415mm	415	4.7	5.3
71446501	14Fr*465mm	465	4.7	5.3
71456501	14Fr*565mm	565	4.7	5.3

## Structure



### Princeton Medical Scientific Co., Ltd

RM503, Building C, No.8 Fukang West Street, Yushu Industrial Park,  
Guangzhou Science City, Luogang District, Guangzhou, 510663, China  
Tel: +86 20 82112462  
Email: pulin@pulinmed.com  
www.pulinmed.com